

Approved Cleaners for Amalgam Separators

- Vacuum lines that discharge amalgam process wastewater to a wastewater facility must not be used with oxidizing cleaners including, but not limited to bleach, chlorine, iodide or peroxide. The pH of the cleaner must be above 6 and below 8. These line cleaners are approved to use with your amalgam separator:

Line Cleaner	Manufacturer
AquaSafe	Compliance Training Partners
Bio-Pure eVacuation System Cleaner	Bio-Pure Products
BioVac	Micrylium
BLAST	Preferred Dist., Inc.
Citrizyme	Pascal Company
DRNA Vac	DRNA
EmPower	Metrex
Enzymax	Hu-Friedy
EZ-Zyme	Miltex
Fresh-Vac	Huntington/Ecolan
Gobble Plus	Germiphene
Green and Clean M2	Metasys
Hayes Evacuation Cleaner	Hayes Hendpiece
Medical Enzyme Detergent	Enzyme Solutions
Neutravac	Biotrol
ProEZ	Cetrol International
Pure-Vac Hg	Sultan Chemists
Sani-Soak Ultra	Enzyme Industries
Sani-TreetPlus	Enzyme Industries
SlugBuster	RAMVAC
SRG Evacuation Line Cleaner	Mydent International
Super Vac 40 Plus	ERP Industries
VacuCleanse	Infection Control Tech
Vacuum Clean	Palermo

- Note City of Escondido makes no recommendations based on product, this list is a general information reference sheet.